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**For women only:** the test should not be performed during or for three days after your menstrual period. The urine sample should not be contaminated with vaginal fluids since this may produce a misleading result.

**Do not make any important medical decision without first referring to your doctor.**

**URIN SAMMELN**

**PROCEDURE**

**ATTENTION:** It's recommended to take a sample of urine for the test in the early morning since it's the most concentrated. The urine used for the test should not come into contact with water from the toilet or any disinfectant or cleaning substances.

**For women only:** the test should not be performed during or for three days after your menstrual period. The urine sample should not be contaminated with vaginal fluids since this may produce a misleading result.

**Do not make any important medical decision without first referring to your doctor.**

**COLLECT URINE:** Collect part of the urine using a clean cup without any residual detergents. Make sure to fill up the cup with urine. - **FIG. A**

**Attention:** The use of urine preservatives is not recommended. If testing cannot be done within an hour after voiding, refrigerate the specimen immediately and let it return to room temperature before testing. Prolonged storage of unpreserved urine at room temperature may result in microbial proliferation with resultant changes in pH. A shift to alkaline pH may cause false positive results with the protein test area. Urine containing glucose may decrease in pH as organisms metabolize the glucose. Contamination of the urine specimen with skin cleansers containing chlorhexidine may affect protein test results.

**PERFORMING THE TEST:**

1) Open the aluminium pouch where indicated and take out the test strip. Do not touch the test fields and dispose of the attached desiccant with your normal household waste. Once opened the pouch, it is recommended to perform the test immediately - **FIG. B**

2) Dip the test strip in the urine sample. **ATTENTION:** Make sure that all four test fields are immersed for about 1-2 seconds. - **FIG. C**

3) Then remove the test strip and wipe off any surplus urine against the rim of the container or with an absorbent material (e.g. a paper towel) to avoid mixing chemicals from adjacent reagent areas.

Place the test strip horizontally on the opened aluminium pouch so as to avoid spilling any urine on any other surfaces. - **FIG. D** **Wait for 2 min.** (do not read results after 3 minutes)

#### RESULT INTERPRETATION


**READ THE RESULT AFTER 2 MINUTES.**

**READ THE RESULT SEPARATELY FOR EACH PARAMETER.**

There are more possible tones and color intensities. Color changes on the edges of the test pads or color changes after more than 3 minutes have to be ignored.

#### NEGATIVE RESULT


The Test field for LEUKOCYTES stayed whitish.  
The Test field for BLOOD stayed mustard yellow.  
The Test field for NITRITE stayed white.  
The Test field for PROTEIN stayed yellowish.

**FIG. E**


**POSITIVE RESULT FOR LEUKOCYTES**

If the color of the test field has changed to purple, then leukocytes have been found in your urine. (see **FIG. E** for the negative result).


**POSITIVE RESULT FOR BLOOD**

If the color of the test field has changed to green (or some green spots appear on the background), then blood has been found in your urine. (see **FIG. E** for the negative result).

**POSITIVE RESULT FOR NITRITE**

If the color of the test field has changed to pink, then nitrites have been found in your urine. (see **FIG. E** for the negative result).

**POSITIVE RESULT FOR PROTEIN**

If the color of the test field has changed to green, then proteins have been found in your urine. (see **FIG. E** for the negative result).

#### PERFORMANCE

**Analytical sensitivity:** Leukocytes: Detects leukocytes as low as 9 white blood cells Leu/ $\mu$ L in clinical urine. Blood: Detects free hemoglobin as low as 0.018 mg/dL or 5 Ery/ $\mu$ L in urine specimens with ascorbic acid content of < 50 mg/dL. Nitrite: Detects sodium nitrite as low as 0.05 mg/dL in urine with a low specific gravity and less than 30 mg/dL ascorbic acid. Protein: Detects albumin as low as 7.5 mg/dL (0.075 g/L).

#### TECHNICAL NOTES ON PARAMETERS

The Test detects LEUKOCYTES, BLOOD, NITRITE and/or PROTEIN in urine.

**LEUKOCYTES:** When taking cephalixin and cephalothin, or high concentration of oxalic acid may also cause test results to be artificially low. Tetracycline may cause decreased reactivity, and high levels of the drug may cause a false negative reaction. High urinary protein may diminish the intensity of the reaction color. The presence of leukocytes in urine is an important symptom of an inflammation of the kidneys and the urinary tract. **BLOOD:** A uniform green color conversion indicates the presence of hemoglobin or hemolyzed erythrocytes; scattered or compacted green spots indicate intact erythrocytes. **NITRITE:** Gram-negative bacteria in urine convert nitrate from food into nitrite. Nitrite reacts with a chemical in the test field and leaves a pink shade. The test result may be distorted if urine does not stay for long in the bladder, due to hunger, a vegetable-free diet or antibiotic treatment and medicines containing phenazopyridine. Comparing the test on a white background may aid in the detection of low nitrite levels, which might otherwise be missed. **PROTEIN:** An indicator on the test field reacts with protein in the urine, changing its colour to green. They may be found where there is inflammation of the bladder or prostate or bleeding in the urinary tract. Infusions containing polyvinylpyrrolidone or medicines containing phenazopyridine may yield a false positive result.





Chemical components in the test fields must be viewed as potentially dangerous substances, although they present no hazard provided that all test components are used in accordance with these instructions.

#### F.A.Q. - QUESTIONS AND ANSWERS

**WHAT SHOULD I DO IF MY TEST RESULT IS POSITIVE?** Remember that a positive result doesn't mean that all four substances have been detected in your urine. Even if your result is positive for just one of them, it is most likely that something is wrong in your urine, even if the reason may not be a urinary infection. Get in touch promptly with your own doctor, who will be able to give a more accurate diagnosis. When you visit your doctor, please take these instructions with you so that he/she will be better informed as to the type of test you have performed.

**WHAT SHOULD I DO IF MY TEST RESULT IS NEGATIVE?** Remember that your test result is only negative if the result on the test field for all four substances is negative. But if you still feel the signs of a UTI or have any other symptom, then contact your own doctor to arrange a more thorough examination.

**WHAT IS THE ACCURACY OF URINARY TRACT INFECTIONS TEST?** One study carried out in parallel with an another test for the urine showed an accuracy of 97.6% for leukocytes, 96% for the blood, 100% for nitrates and 88% for proteins.

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#### F.A.Q. - QUESTIONS AND ANSWERS

**WHAT SHOULD I DO IF MY TEST RESULT IS POSITIVE?** Remember that a positive result doesn't mean that all four substances have been detected in your urine. Even if your result is positive for just one of them, it is most likely that something is wrong in your urine, even if the reason may not be a urinary infection. Get in touch promptly with your own doctor, who will be able to give a