









READING THE RESULTS

	Normal Two colored lines appear. Both T (Test) and C (Control) line appear. This result means that the Ferritin concentration in blood is normal.
	Abnormal One colored line appears. Only control line (C) appears. This result means that the ferritin concentration in blood is too low. You should consult physician.
	Invalid Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor

LIMITATIONS

- The Iron Fer Test provides only a qualitative analytical result. A secondary analytical method must be used to obtain a confirmed result.
- It is possible that technical or procedural errors, as well as other interfering substances in the whole blood specimen may cause erroneous results.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- Other clinically available tests are required if questionable results are obtained.

PERFORMANCE CHARACTERISTICS

The specimen correlation used a specimen number (n) equal to 102 specimens, including 79 normal whole blood specimens and 23 abnormal whole blood specimens were confirmed by CLIA. The result demonstrated showed that the abnormal coincidence rate is 91.3%, the normal coincidence rate is 96.2% and the total coincidence rate is 95.1%.
Iron Fer Test Result

Method	CLIA			Total Results
	Results	Abnormal	Normal	
Iron Fer Test	Abnormal	21	3	24
	Normal	2	76	78
Total Results		23	79	102

Abnormal coincidence rate= $\frac{21}{(21+2)} \times 100\% = 91.3\%$ Normal coincidence rate= $\frac{76}{(3+76)} \times 100\% = 96.2\%$
Total coincidence rate= $\frac{(21+76)}{(21+3+2+76)} \times 100\% = 95.1\%$

Accuracy

The Iron Fer Test has been compared with a leading commercial Ferritin CLIA test. The correlation between these two systems is over 95.0%.

Intra-Assay

Within-run precision has been determined by using 10 replicates of three specimens: 0 ng/mL, 30 ng/mL and 100 ng/mL specimens. The specimens were correctly identified > 99% of the time.

Inter-Assay

An evaluation between-run precision has been determined by 10 independent assays on the same 3 specimens: 0 ng/mL ferritin, 30 ng/mL ferritin, 100 ng/mL ferritin standard sample. Three different lots of the Iron Fer Test have been tested using these specimens. The specimens were correctly identified >99% of the time.

Analytical Sensitivity

The Iron Fer Test can detect levels of ferritin in human fingerstick blood as low as 30 ng/mL.

Cross-Reactivity

An evaluation was performed to determine the cross reactivity and interferences of Iron Fer Test. There is no cross reactivity with HAMA, RF, Human serum albumin, human AFP, Ferric Chloride, human transferrin and human hemoglobin.

EXTRA INFORMATION

1. How does the ferritin test work?

Ferritin is a protein and the primary form of iron stored inside cells. An abnormal result means that the ferritin concentration in blood is lower than 30 ng/mL and a possible iron deficiency.

2. When should the test be used?

The Iron Fer Test can be performed in case of symptoms like paleness, feeling tired, headaches, faster heartbeat or shortness of breath during exercise; mainly, if woman, when pregnant or in case of excessive bleeding during periods. The test can be performed anytime of the day, but must not be performed in case of disease, acute inflammations or in case of spleen or liver injury. Abnormal results can be obtained even in case of no iron deficiency situation.

3. Can the result be incorrect?

The results are accurate as long as the instructions are carefully respected. Nevertheless, the result can be incorrect if the ferritin test gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient. The capillary dropper provided in the box allows making sure the collected blood volume is correct. Besides, due to immunological principles involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.

4. What is the line that appears under the C (control) line?

When this line appears, it only means that the test is performing well.

5. If i read the result after 10 minutes, will the result be reliable?

No. The result should be read at 5 minutes after adding the buffer. The result is not reliable after 10 minutes.

6. What do I have to do if the result is abnormal?

If the result is abnormal, it means that the ferritin level is lower than the normal (30 ng/mL) and that you should consult the physician and show the test result to him/her. Then, the physician will decide whether additional analysis should be performed.

7. What do I have to do if the result is normal?

If the result is normal, it means that the ferritin level is higher than 30 ng/mL and is within the normal range. However, if the symptoms persist, it is recommended to consult a physician.

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